



LIQUID LAUNDRY BLOOD REMOVER

PRODUCT USES

Blood stain remover used to pre-treat blood stained fabrics prior to washing.

OPERATOR BENEFITS

- Reduces purchases of replacement linen.
- Quickly penetrates tough to remove blood stains.

DIRECTIONS FOR USE

1. Apply full-strength Blood Remover directly on the blood stains.
2. Gently rub the stained areas to ensure complete penetration of the stains.
3. Wash normally using the appropriate laundry detergent and wash formula.

PRODUCT SPECS

Pack sizes:	6/1 Qt.	Color:	Clear straw
Characteristics:	Non-viscous liquid	Odor:	Slight amonia
Rinsability:	Excellent	Foam:	None
Flammability:	NA	pH in conc. form:	7.5 – 8.5
		Phosphorous content:	0.75%

IMPORTANT OPERATOR INFORMATION

MSDS:

Refer to the Material Safety Data Sheet (MSDS) supplied by your distributor for information on emergency response, safe usage, hazards, first aid and proper product handling.

Emergency Response Phone:

1-800-424-9300

Statement of

Practical Treatment:

If swallowed: Dilute by giving large amounts of water and call a physician.

For eye or skin contact: Flush with cool water for 15 minutes and call a physician.

Product Precautions:

CAUTION. Harmful if swallowed. Injurious to eyes and skin.

DANGER: KEEP OUT OF REACH OF CHILDREN.

Shipping Classification:

Not Regulated

