



STAINWORX



Blood Buster

PRODUCT USES	Blood stain remover used to pre-treat blood stained fabrics prior to washing.	
OPERATOR BENEFITS	<ul style="list-style-type: none">• Reduces purchases of replacement linen.• Quickly penetrates tough to remove blood stains.	
DIRECTIONS FOR USE	<ol style="list-style-type: none">1. Apply full-strength Blood Remover directly on the blood stains.2. Gently rub the stained areas to ensure complete penetration of the stains.3. Wash normally using the appropriate laundry detergent and wash formula.	
PRODUCT SPECS	Pack size: 6/1 QT Characteristics: Non-viscous liquid Rinsability: Excellent Flammability: NA	Color: Clear straw Odor: Slight amonia Foam: None pH in conc. form: 7.5 – 8.5 Phosphorous content: 0.75%

IMPORTANT OPERATOR INFORMATION

SDS:	Refer to the Safety Data Sheet (SDS) supplied by your distributor for information on emergency response, safe usage, hazards, first aid and proper product handling.
Emergency Response Phone:	1-800-424-9300
Statement of Practical Treatment:	If swallowed: Dilute by giving large amounts of water and call a physician. For eye or skin contact: Flush with cool water for 15 minutes and call a physician.
Product Precautions:	CAUTION. Harmful if swallowed. Injurious to eyes and skin. DANGER: KEEP OUT OF REACH OF CHILDREN.
Shipping Classification:	Not Regulated